



FLOW CYTOMETRY TROUBLESHOOTING GUIDE



Quickly diagnose and fix common experimental issues

WHAT PROBLEM ARE YOU SEEING?

"Start here first"

"Most common issue"

"Often overlooked"

WEAK OR NO SIGNAL

CAUSES

- Low antigen expression
- Antibody too dilute
- Low detector sensitivity

FIXES

- Increase antibody concentration
- Use brighter fluorophore
- Adjust PMT voltage

HIGH BACKGROUND

CAUSES

- Non-specific binding
- Excess antibody
- Dead cells present

FIXES

- Add Fc blocking step
- Reduce antibody concentration
- Use viability dye

POOR RESOLUTION

CAUSES

- Incorrect compensation
- Spectral overlap
- Suboptimal voltages

FIXES

- Use single-color controls
- Adjust compensation matrix
- Optimize panel design

DEBRIS / CLOGGING

CAUSES

- Cell aggregates
- Unfiltered samples
- Instrument blockage

FIXES

- Filter samples before run
- Wash cells thoroughly
- Perform instrument cleaning

UNEXPECTED POPULATIONS

CAUSES

- Improper gating
- Doublets or aggregates
- Sample contamination

FIXES

- Apply standard gating strategy
- Exclude doublets
- Verify controls and labeling

BEST PRACTICES

Prevent issues before they start



TITRATE all antibodies

INCLUDE proper controls

MAINTAIN instrument QC

FILTER all samples

PROTECT from light

