


APC Streptavidin

Catalog Number:	306301, 306302
Size:	25 ug, 100 ug
Target Name:	Streptavidin
Regulatory Status:	RUO

PRODUCT DETAILS

Application:	Flow Cytometry
Format:	PE
Formulation:	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide and 0.2% (w/v) BSA
Protein Concentration:	0.2 mg/mL
Recommended Usage:	For flow cytometric staining, it is recommended to use less than 0.1 ug of this reagent per 0.5-1.0 million cells in a 100 μ L volume. Optimal reagent performance should be determined by titration for each specific application. PE has an excitation max at 565 nm and an emission max at 575 nm.
Excitation Laser:	Red Laser (633 nm)

PRODUCT DATA



A-431 cells stained either Biotin Anti-Human EGFR clone 528 (color-filled histogram) or an isotype control (gray histogram), followed by APC Streptavidin.