

InnoCyto Inc.

15375 Barranca Pkwy, Suite I-103 Irvine, CA 92618

Technical Data Sheet

Anti-Human CD34

Catalog Number: 107001

Size: 100 ug

Target Name: CD34, Gp105-120, My10

Regulatory Status: RUO

Product Details

Clone: 034AB
Application: FC
Reactivity: Human
Format: Purified
Isotype: Rabbit IgG

Antibody Type: Monoclonal

Formulation: Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide

Protein Concentration: 0.5 mg/mL

Storage&Handling: The antibody solution should be stored between 2°C and 8°C

Recommended Usage: For flow cytometric staining, it is recommended to use less than 0.2 ug of this reagent per 0.5-1.0 million cells in a 100 μ L volume. Optimal reagent performance should be

determined by titration for each specific application

Release Date: Dec-25

Background Information


CD34 (gp105–120) is a type I monomeric sialomucin-like glycophosphoprotein with a molecular weight of 105–120 kDa, expressed on hematopoietic stem and progenitor cells, bone marrow stromal cells, capillary endothelial cells, embryonic fibroblasts, and some neural tissues. CD34 serves as a widely used stem cell marker and mediates cell adhesion and lymphocyte homing through interactions with L-selectin and E-selectin ligands. CD34 expression is dynamic and reversible, reflecting specific states of hematopoietic development, including proliferation, differentiation, and altered adhesion properties. Multiple epitopes have been identified based on differential enzymatic sensitivity. Clinically, CD34 is expressed in ~40% of acute myeloid leukemias and 65% of pre-B acute lymphoblastic leukemias, but only 1–5% of acute T-lymphoid leukemias. It is also used to assess stem cell counts for transplantation and can provide prognostic information in certain cancers.

Product Data



InnoCyto Inc.

15375 Barranca Pkwy, Suite I-103 Irvine, CA 92618



Human peripheral blood leukocytes stained with iF647 Anti-human CD45 and either purified Anti-Human CD34 clone 034AB (right panel) or an isotype control (left panel), followed by PE anti-Ratbbit IgG 2nd antibody.